

Pocket design with simple operation
Daily self-check at home



Wireless connection to smartphone via APP
Real-time display and data storage



Cloud server for data storage
Records within reach anywhere



Data to be remotely reviewed &
analyzed by the medical team



PC-80B

Easy ECG Monitor

 CREATIVE MEDICAL  CE ₀₁₂₃ FDA 510(k)

Shenzhen Creative Industry Co., Ltd.


















Floor 5, BLD 9, Baiwangxin High-Tech Industrial Park, Songbai Road, Xili Street
Nanshan District, 518110 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Tel: +86 755 2643 3514 Fax: +86 755 2643 0390

www.creative-sz.com info@creative-sz.com

 CREATIVE MEDICAL

17 types of measurement


No irregularity found 	Suspected occasional short beat interval 	Suspected fast beat with baseline wander 
Suspected a little fast beat 	Suspected irregular beat interval 	Suspected slow beat with baseline wander 
Suspected fast beat 	Suspected fast beat with short beat interval 	Suspected occasional short beat interval with baseline wander 
Suspected short run of fast beat 	Suspected slow beat with short beat interval 	Suspected irregular beat interval with baseline wander 
Suspected a little slow beat 	Suspected slow beat with irregular beat interval 	Poor Signal, measure again 
Suspected slow beat 	Waveform baseline wander 	

Remark: The above measurement results are for user's reference only



Features:


- Palm design for convenient daily measurement at home
- Palm/Leg/Chest three quick measurement methods within 30 seconds
- 22-hour storage of continuous measurement
- 17 types of preset ECG waveform analysis
- Android/IOS application for real-time display and storage
- Data management software for data review and analysis
- Communication protocol, SDK available for integration




ECG Viewer Manager

Technical Specifications


Number of channel	Single channel
Electrodes	3 embedded metal electrodes or external ECG leadwire
Display	2.8 inch color LCD screen
HR measuring range	30bpm~240bpm
HR measuring accuracy	±2bpm or ±2%
Display scale	5.0mm/mV ±10%
Sweeping speed	20mm/s ±10%
Measurement time	30s quick measurement or continuous measurement
Measurement result	17 types
Data storage	2700 groups record or 22 hours trend
Data upload	By USB with PC analysis software
Data transmission	By wireless real-time transfer
Dimension	125mm(L) x 70mm(W) x 21.5mm(H)
Net weight	106g (including two AAA batteries)
Power supply	DC 3V(2xAAA battery)




Palm measurement



Leg measurement



Chest measurement



Lead Wire measurement

