

# LMT 10 Digital Thermometer

LMT 10 digital thermometer has a pure white light body. The light weight and small size make the product easy to use and carry. The digital thermometer can measure the human body temperature quickly and accurately with the buzzing reminder. In addition to high measurement accuracy, LMT 10 is also a device with IP22 waterproof level.

Lightweight

Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard

**Division Value:**

0.1°C / 0.1°F

**Measurement Range:**

32.0 °C to 42.9 °C (89.6°F to 109.2°F)

**Product Size:**

130.5(length) mm X 20.5  
(width) mm X 11.4(height) mm

**Accuracy:**

(35°C ~ 42.0°C) ±0.1°C;  
(32.0 ~ 34.9 °C and 42.1 ~ 42.9 °C) ± 0.2 °C

**Buzzing Reminder:**

Yes

**Product Weight:**

12g (Including Battery)

**Can Be Switched  
Between °C and °F:**

Yes

**Water Proof:**

Yes

**Battery Lifetime:**

Approx. 10000 Measurements  
(Using a New Battery)

**Battery:**

LR41 (1.5V) / SR41 (1.55V)

**Operating Conditions:**

5-40 °C

**Automatic Shutdown Time:**

8 Minutes 40 Seconds ± 1 Minute

**Storage / Transport Conditions:**

-20 ~ +55 °C

**Normal Body Temperature Range:**

Armpit temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C;  
Oral Temperature: 35.9°C-36°C



**Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.**

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278   Fax: +86 0755-86952278   ✉ E-mail: marketing@lepu-medical.com

# LMT 11A Digital Thermometer

LMT 11A digital thermometer has the light body with blue and white. Its light weight makes it easy to use and carry. The intelligent reminder function makes it convenient for people to use it, and LMT 11A can be also switched the temperature between Fahrenheit and Celsius.

Lightweight















Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard

 <b>Division Value:</b> 0.1°C /0.1°F	 <b>Measurement Range:</b> 32.0 °C to 42.9 °C (89.6°F to 109.2°F)
 <b>Product Size:</b> 131.8(Length) mm X 21.1 (Width) mm X 11.9(Height) mm	 <b>Accuracy:</b> (35°C ~ 42.0°C) ±0.1°C; (32.0 ~ 34.9 °C and 42.1 ~ 42.9 °C) ± 0.2 °C
 <b>Buzzing Reminder:</b> Yes	 <b>Product Weight:</b> 12g (Including Battery)
 <b>Can Be Switched Between °C and °F:</b> Yes	 <b>Water Proof:</b> Yes
 <b>Battery Lifetime:</b> Approx. 10000 Measurements (Using a New Battery)	 <b>Battery:</b> LR41 (1.5V) / SR41 (1.55V)
 <b>Operating Conditions:</b> 5-40 °C	 <b>Automatic Shutdown Time:</b> 8 Minutes 40 Seconds ± 1 Minute
 <b>Storage / Transport Conditions:</b> -20 ~ +55 °C	
 <b>Normal Body Temperature Range:</b> Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C; Oral Temperature: 35.9°C-36°C	



**Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.**

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278    📠 Fax: +86 0755-86952278    ✉ E-mail: marketing@lepu-medical.com

# LMT 11B Digital Thermometer

LMT 11B digital thermometer has the light body with blue and white. Its temperature measuring probe is soft, so the LMT 11B is also suitable for the elderly and babies. The intelligent reminder function makes it convenient for people to use it, and LMT 11B can be also switched the temperature between Fahrenheit and Celsius.

**Lightweight**














**Reminder**

**Fast & Accurate**



# Performance Parameters

Standard Type: Medical Standard

 <b>Division Value:</b> 0.1°C / 0.1°F	 <b>Measurement Range:</b> 32.0 °C to 42.9 °C (89.6°F to 109.2°F)
 <b>Product Size:</b> 131.8(Length) mm X 21.1 (Width) mm X 11.9(Height) mm	 <b>Accuracy:</b> (35°C ~ 42.0°C) ±0.1°C; (32.0 ~ 34.9 °C and 42.1 ~ 42.9 °C) ±0.2 °C
 <b>Buzzing Reminder:</b> Yes	 <b>Product Weight:</b> 12g (Including Battery)
 <b>Can Be Switched Between °C and °F:</b> Yes	 <b>Water Proof:</b> Yes
 <b>Battery Lifetime:</b> Approx. 10000 Measurements (Using a New Battery)	 <b>Battery:</b> LR41 (1.5V) / SR41 (1.55V)
 <b>Operating Conditions:</b> 5-40 °C	 <b>Automatic Shutdown Time:</b> 8 Minutes 40 Seconds ± 1 Minute
 <b>Storage / Transport Conditions:</b> -20 ~ +55 °C	
 <b>Normal Body Temperature Range:</b> Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C; Oral Temperature: 35.9°C-36°C	



**Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.**

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278    📠 Fax: +86 0755-86952278    ✉ E-mail: marketing@lepu-medical.com



# LMT 11C Digital Thermometer

LMT 11C digital thermometer has a lightweight body with gray and white. Buzzing reminder and water-proof function are its two highlights. Whether at home or in the hospital, people can measure body temperature quickly and accurately.

Lightweight















Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard

 <b>Division Value:</b> 0.1°C / 0.1°F	 <b>Measurement Range:</b> 32.0 °C to 42.9 °C (89.6°F to 109.2°F)
 <b>Product Size:</b> 131.8(length) mm X 21.1 (width) mm X 11.9(height) mm	 <b>Accuracy:</b> (35°C ~ 42.0°C) ±0.1°C; (32.0 ~ 34.9 °C and 42.1 ~ 42.9 °C) ± 0.2 °C
 <b>Buzzing Reminder:</b> Yes	 <b>Product Weight:</b> 12g (Including Battery)
 <b>Can Be Switched Between °C and °F:</b> Yes	 <b>Water Proof:</b> Yes
 <b>Battery Lifetime:</b> Approx. 10000 Measurements (Using a New Battery)	 <b>Battery:</b> LR41 (1.5V) / SR41 (1.55V)
 <b>Operating Conditions:</b> 5-40 °C	 <b>Automatic Shutdown Time:</b> 8 Minutes 40 Seconds ± 1 Minute
 <b>Storage / Transport Conditions:</b> -20 ~ +55 °C	
 <b>Normal Body Temperature Range:</b> Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C; Oral Temperature: 35.9°C-36°C	



**Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.**

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278    📠 Fax: +86 0755-86952278    ✉ E-mail: [marketing@lepu-medical.com](mailto:marketing@lepu-medical.com)




# LMT 12A Digital Thermometer

LMT 12A digital thermometer can measure the human body temperature quickly and accurately in 20 seconds. It has the characteristics of convenient reading, short measurement time, high measurement accuracy, buzzer reminder, and no mercury. It is especially suitable for use in hospitals, homes and other environments.

Lightweight
















Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard

 <b>Division Value:</b> 0.1°C /0.1°F	 <b>Measurement Range:</b> 32.0 °C to 42.9 °C (89.6°F to 109.2°F)
 <b>Product Size:</b> 133.5(length) mm X 32.0 (width) mm X 15.0(height) mm	 <b>Accuracy:</b> (35°C ~ 42.0°C) ±0.1°C; (32.0 ~ 34.9 °C and 42.1 ~ 42.9 °C) ± 0.2 °C
 <b>Buzzing Reminder:</b> Yes	 <b>Approx. Measuring Time:</b> 20 Seconds
 <b>Battery Lifetime:</b> Approx. 10000 Measurements (Using a New Battery)	 <b>Product Weight:</b> 22g (Including Battery)
 <b>Can Be Switched Between °C and °F:</b> Yes	 <b>Water Proof:</b> Yes
 <b>Operating Conditions:</b> 5-40 °C	 <b>Battery:</b> CR2032 (3.0V)
 <b>Storage / Transport Conditions:</b> -20 ~ +55 °C	 <b>Automatic Shutdown:</b> Yes
 <b>Normal Body Temperature Range:</b> Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C; Oral Temperature: 35.9°C-36°C	



**Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.**

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278    📠 Fax: +86 0755-86952278    ✉ E-mail: [marketing@lepu-medical.com](mailto:marketing@lepu-medical.com)

# LMT 12B Digital Thermometer

LMT 12B digital thermometer can measure the human body temperature quickly and accurately in 20 seconds. It has the soft temperature probe, the large screen, the buzzing reminder function, high measurement accuracy, and no mercury. It's especially suitable for use in hospitals, homes and other environments.

Lightweight
















Reminder

Fast & Accurate



# Performance Parameters

Standard Type: Medical Standard

 <b>Division Value:</b> 0.1°C / 0.1°F	 <b>Measurement Range:</b> 32.0 °C to 42.9 °C (89.6°F to 109.2°F)
 <b>Product Size:</b> 133.5(length) mm X 32.0 (width) mm X 15.0(height) mm	 <b>Accuracy:</b> (35°C ~ 42.0°C) ±0.1°C; (32.0 ~ 34.9 °C and 42.1 ~ 42.9 °C) ±0.2 °C
 <b>Buzzing Reminder:</b> Yes	 <b>Approx. Measuring Time:</b> 20 Seconds
 <b>Battery Lifetime:</b> Approx. 10000 Measurements (Using a New Battery)	 <b>Product Weight:</b> 22g (Including Battery)
 <b>Can Be Switched Between °C and °F:</b> Yes	 <b>Water Proof:</b> Yes
 <b>Operating Conditions:</b> 5-40 °C	 <b>Battery:</b> CR2032 (3.0V)
 <b>Storage / Transport Conditions:</b> -20 ~ +55 °C	 <b>Automatic Shutdown:</b> Yes
 <b>Normal Body Temperature Range:</b> Armpit Temperature: 36.0°C-37.3°C; Rectal Temperature: 36.5°C-37.6°C; Oral Temperature: 35.9°C-36°C	



**Shenzhen LEPU Intelligent Medical Equipment Co., Ltd.**

📍 Add : North Side of Floor 3, BLD 9 BaiWangxin High-Tech Industrial Park, Songbai Road, Xili Street, Nanshan District 518055 Shenzhen, Guangdong, CHINA

☎ Tel: +86 0755-86952278    📠 Fax: +86 0755-86952278    ✉ E-mail: [marketing@lepu-medical.com](mailto:marketing@lepu-medical.com)